


# Keller Meditech Installation

Log into the Sullivan VPN from our webpage <http://remote.helenkeller.com>



The screenshot shows the Helen Keller Hospital Remote Access Web Page. At the top left is the Helen Keller Hospital logo. The main heading is "Helen Keller Hospital Remote Access Web Page". Below this is a grid of six links:

<b><a href="#">SULLIVAN VPN</a></b> Establish VPN Connection <a href="#">Instruction on establishing a VPN Connection</a> <a href="#">Instruction on Installing Meditech</a>	<b><a href="#">WEB PORTAL</a></b> Access Keller's Web Based Applications <a href="#">PC Requirements for the Keller Web Portal</a>	<b><a href="#">EHR PORTAL</a></b> Physician Electronic Health Records Esign <a href="#">Instructions for Remote Access and EHR Portal</a>
<b><a href="#">OWA EMAIL</a></b> Outlook Webmail Access <a href="#">Outlook Web Access Quick Start Guide</a>	<b><a href="#">INFORMATION SYSTEMS</a></b> Keller Information Systems Department	

A red arrow points to the "SULLIVAN VPN" link.

Login with your VPN credentials, click on the **HOME** button on the right hand column. Then click on the **Meditech Install Program**.




The screenshot shows the Helen Keller Sullivan Login page. At the top left is the Helen Keller Hospital logo. The main heading is "Helen Keller Sullivan Login". Below this is a navigation bar with "Home" and "AnyConnect" buttons. A red arrow points to the "Home" button. Below the navigation bar is a search bar labeled "Address". Below the search bar is a photograph of the hospital building. Below the photograph is the heading "Sullivan AnyConnect VPN Web Portal". Below this heading is a table with three rows:

<a href="#">GE Centricity Pictorial Viewer Patch</a>	Launch this patch if you are having problems viewing images in GE Centricity
<a href="#">Meditech Installation Instructions</a>	Instructions on how to install and setup Meditech
<a href="#">Meditech Install Program</a>	Setup application for Meditech

A red arrow points to the "Meditech Install Program" link.


A short installation box will display followed by a prompt that states that the application is up to date.




Double click on the icon **MAGIC Workstation** icon on your desktop.




A connection will have to be defined the first time MAGIC Workstation (Meditech) is launched. To create a connection click ADD.




The connection name is **KELLER**, the connection Type is **Telnet**. Click **OK**.



Setup the KELLER connection using **MAGIC** terminal emulation with a Host Name/IP Address **128.1.0.1** and Telnet Port of **23**. When you are finished hit **OK**.



Click **CONNECT** on the following screen to start your connection to Meditech.



**Please make sure you disconnect your VPN Session when you are done using Meditech.**